



Introduction to PIV

Part 1
Collecting Good PIV Data
TSI Webinar Training Series

March 3, 2010
Matt Stegmeir
Tyson Strand



Topics Covered



- **PIV Fundamentals (2D PIV)**
 - Cross Correlation
 - Rules of Thumb
- **Experimental Considerations**
- **Measurement Setup**
 - Calibration
 - Seed Density
 - Optimizing Δt
- **Introduction to Insight 3G**

PIV Fundamentals





- **PIV measures the displacement of tracer seed particles in a flowfield to determine the 2D velocity field**
 - **At time t_1**
 - **Pulsed laser sheet illuminates a planar region of the flow**
 - **Particles are imaged on the camera (Frame A)**
 - **At time $t_1 + \Delta t$**
 - **A second image (Frame B) is taken of a second light sheet**
 - **Statistical (Cross-Correlation) methods are used to determine the particle displacement over the time Δt , and thus the local velocity**


Cross Correlation




- All particles look alike, so it is hard to find the 'same' particle in both Frame A and B (PTV)
- Instead, PIV uses a statistical approach to find the most likely displacement of a group of particles
- Frame A is broken up into a grid of 'interrogation regions'
- The group of particles in the interrogation region creates a fairly unique 'fingerprint' that we can look for in both frames (PIV)



Cross Correlation




Cross Correlation





Original position
and particle images
from Frame A

Cross Correlation




Particle Images in
Frame B

Cross Correlation




A 'Spot Mask' (above) is scanned through the search area in Frame B, and a correlation value is calculated at each position

Cross Correlation




The value of the correlation function will be maximum when the 'fingerprint' in the Spot Mask is identified in Frame B

Cross Correlation




The X and Y displacements of the particles are determined by the offset of the interrogation regions

Cross Correlation



Since we know the time separation (Δt) between the 2 laser pulses very accurately, we can measure the group velocity as the displacement / Δt , and assign a single velocity vector

Cross Correlation



The process is repeated for each interrogation region in Frame A, resulting in a 2-Dimensional velocity field for the imaged region



PIV Rules of Thumb




- The raw data of PIV measurements are particle images
- The key to good measurements is good raw data
- In the Ideal Case
 - Seed particle images are 3 – 5 pixels in diameter
 - 5 – 15 particles per interrogation region
 - Maximum particle displacement approximately 25% of size of interrogation region
- It is not always possible to satisfy these in all measurements, but they are good experimental goals

Experiment Design



- Important considerations:
 - What size field of view?
 - What desired spatial resolution?
 - Are the above 2 realistic together?
 - Appropriate seeding
 - Optical access
 - Camera can view measurement region at 90 degrees
 - View is not distorted
 - What does it look like when you view the measurement region from the perspective of the camera?



With a 4MP (2k x 2k) Camera

- 100 microns per pixel
- Spatial resolution
 - 3.2 mm at 32 pixel IR
 - 6.4 mm at 64 pixel IR

Experimental Setup



- Assuming laser is set up to illuminate desired plane in flow, camera at 90 degrees to light sheet with good optical access (no distortions)
- Initial setup steps should be performed in 'free, continuous' mode in the Insight 3G software
 - Using room light, no laser
- For final setup, the best practice is to focus the camera on seed particles in the flow
 - This assures us that our focus is optimized to the laser sheet



Experimental Hints



- **When setting up the cameras to view laser light scattered from particles**
 - **Start on a large f-number (small aperture) and fairly low laser energy**
 - **Increase alternately until seed particles are well-illuminated with minimal pixel saturation**
 - **WARNING: IF LARGE REGIONS OF SATURATED (PINK) PIXELS APPEAR, STOP CAPTURE IMMEDIATELY**
 - **Laser light can damage camera pixels**
 - **Reduce laser energy and/or increase f-number**

Calibration




- **What are we ‘calibrating’?**
 - We need to show the software how to convert from pixel units on the images, to the physical units in our flow space
- **We are considering the simplest case in 2D PIV, where we view the light sheet at 90 degrees with good optical access**
- **Off-axis viewing and optical distortions are sometimes unavoidable, and appropriate corrections can be made**
 - In these cases, we need to give the software more information to be able to convert pixels to mm over the entire image region
 - Addressed in future webinars

Calibration Procedure



- Assumption: Cameras focused on laser light sheet
- Calibration Steps
 1. Capture and save image of “Ruler” in plane of laser light sheet
 2. Create a calibration file with the saved image
 3. Measure across a known distance in the calibration image




Seed Density



- The true test of appropriate seed particle density is examining the Particle Images


- Field of View is known
- Desired Spatial Resolution known
- Do I have enough particles in each interrogation region?



Timing Setup



- The Timing Setup Window in Insight 3G



Optimizing Δt




- The next step to Good PIV Data is setting an appropriate value for Δt
- Δt determines the particle image displacement in PIV images
- We must select appropriate Δt so that particle image displacement is consistent with experimental objectives
 - Field of view
 - Spatial resolution
 - Appropriate limits on particle displacement

Optimizing Δt




- The first step in optimizing Δt is developing “the eye”
- Can you see the displacement?
- Although qualitative, this is a critical step in optimizing a PIV measurement
 - If displacement appears random, reduce Δt
 - If there is little / no displacement, increase Δt



Optimizing Δt




- There are several ways to 'measure' the displacement for a more quantitative assessment
 - Zoom into individual particles
 - Perform 'Point Processing' and assess displacement of individual spots
- These will be illustrated in the following demonstration of the Insight 3G Software





Introduction to INSIGHT 3G

- Tricks and techniques for helping you to optimize your measurements using Insight 3G
 - Calibration
 - Timing Setup
 - Assessing Δt
 - Image Zoom
 - Point Process



Thank you!



- Questions?

www.tsi.com/fluidmechanics

fluid@tsi.com